

ARG65848 anti-Dendritic Cells antibody [MRC OX-62] (FITC)

Package: 100 µl
Store at: 4°C

Summary

Product Description	FITC-conjugated Mouse Monoclonal antibody [MRC OX-62] recognizes Dendritic Cells
Tested Reactivity	Rat
Tested Application	FACS
Specificity	This antibody only labeled cells characteristic of veiled cells; the antigen recognized by this antibody is expressed at low to moderate levels by a subset of spleen, thymus and bone marrow dendritic cells.
Host	Mouse
Clonality	Monoclonal
Clone	MRC OX-62
Isotype	IgG1
Target Name	Dendritic Cells
Species	Rat
Immunogen	Rat veiled (dendritic) cells obtained from lymph.
Conjugation	FITC

Application Instructions

Application table	Application	Dilution
	FACS	Assay-dependent
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

Properties

Form	Liquid
Purification	Purified.
Buffer	PBS, 0.02% Sodium azide and BSA.
Preservative	0.02% Sodium azide
Stabilizer	BSA
Storage instruction	Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.