

ARG66312 Goat anti-Rabbit IgG antibody (TAMRA)

Package: 500 µl
Store at: 4°C

Summary

Product Description	TAMRA-conjugated Goat Polyclonal antibody recognizes Rabbit IgG
Tested Reactivity	Rb
Tested Application	FACS, ICC/IF
Host	Goat
Clonality	Polyclonal
Isotype	IgG
Target Name	IgG
Species	Rabbit
Immunogen	Rabbit IgG
Conjugation	TAMRA

Application Instructions

Application table	Application	Dilution
	FACS	1:400 - 1:1500
	ICC/IF	1:400 - 1:1500
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

Properties

Form	Liquid
Purification	Purified.
Buffer	PBS (pH 7.4), 0.01% Thimerosal and 1% BSA.
Preservative	0.01% Thimerosal
Stabilizer	1% BSA
Concentration	1 mg/ml
Storage instruction	Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.