

ARG66354 anti-Cytokeratin 10 antibody [SQab18104]

Package: 100 µl
Store at: -20°C

Summary

Product Description	Recombinant Rabbit Monoclonal antibody [SQab18104] recognizes Cytokeratin 10
Tested Reactivity	Hu
Tested Application	IHC-P
Host	Rabbit
Clonality	Monoclonal
Clone	SQab18104
Isotype	IgG
Target Name	Cytokeratin 10
Species	Human
Immunogen	Synthetic peptide within aa. 484-584 of Human Cytokeratin 10.
Conjugation	Un-conjugated
Alternate Names	KPP; K10; CK-10; BIE; Keratin, type I cytoskeletal 10; Cytokeratin-10; CK10; BCIE; Keratin-10; EHK

Application Instructions

Application table	Application	Dilution
	IHC-P	1:100 - 1:200
Application Note	IHC-P: Antigen Retrieval: Heat mediated was performed using Tris/EDTA buffer (pH 9.0), primary antibody incubate at RT for 30 min. * The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	


Properties

Form	Liquid
Purification	Purification with Protein A.
Buffer	PBS, 0.01% Sodium azide, 40% Glycerol and 0.05% BSA.
Preservative	0.01% Sodium azide
Stabilizer	40% Glycerol and 0.05% BSA
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Gene Symbol	KRT10
Gene Full Name	keratin 10, type I
Background	<p>Cytokeratin 10 is a member of the type I (acidic) cytokeratin family, which belongs to the superfamily of intermediate filament (IF) proteins. Keratins are heteropolymeric structural proteins which form the intermediate filament. These filaments, along with actin microfilaments and microtubules, compose the cytoskeleton of epithelial cells. Mutations in this gene are associated with epidermolytic hyperkeratosis. This gene is located within a cluster of keratin family members on chromosome 17q21. [provided by RefSeq, Jul 2008]</p>
Function	<p>Cytokeratin 10 plays a role in the establishment of the epidermal barrier on plantar skin.</p> <p>(Microbial infection) Acts as a mediator of S.aureus adherence to desquamated nasal epithelial cells via clfB, and hence may play a role in nasal colonization.</p> <p>(Microbial infection) Binds S.pneumoniae PsrP, mediating adherence of the bacteria to lung cell lines. Reduction of levels of KRT10 keratin decrease adherence, overexpression increases adherence. Neither protein has to be glycosylated for the interaction to occur. [UniProt]</p>
Highlight	<p>Related products: Anti-Rabbit IgG secondary antibodies;</p> <p>Related news: Cancer Pathology Markers (SQ clones)</p>
Calculated Mw	59 kDa

Images



ARG66354 anti-Cytokeratin 10 antibody [SQab18104] IHC-P image

Immunohistochemistry: Formalin-fixed and paraffin-embedded cervix tissue stained with ARG66354 anti-Cytokeratin 10 antibody [SQab18104]. Antigen Retrieval: Heat mediated was performed using Tris/EDTA buffer (pH 9.0).