

## ARG67148 Mouse IgG2b Kappa Isotype Control antibody

Package: 1 ml  
Store at: -20°C

### Summary

Product Description	Mouse Monoclonal antibody as a negative control antibody for Mouse IgG2b Kappa. ARG67148 Mouse IgG2b Kappa Isotype Control antibody is Glycerol-free and Carrier-free.
Tested Application	ELISA, FACS, ID, IE, IP, WB
Host	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Target Name	Mouse IgG2b Kappa
Conjugation	Un-conjugated

### Application Instructions

Application table	Application	Dilution
	ELISA	Assay-dependent
	FACS	Assay-dependent
	ID	Assay-dependent
	IE	Assay-dependent
	IP	Assay-dependent
	WB	Assay-dependent

**Application Note** \* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.

### Properties

Form	Liquid
Buffer	PBS (pH 7.4) and 0.03% Proclin-300
Preservative	Proclin-300
Concentration	1 mg/ml
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.