

ARG67203 anti-Canine Coronavirus N Protein antibody [2B8]

Package: 100 µg
Store at: -20°C

Summary


Product Description	Mouse Monoclonal antibody [2B8] recognizes Canine Coronavirus N Protein
Tested Reactivity	Virus
Tested Application	ICC/IF, WB
Host	Mouse
Clonality	Monoclonal
Clone	2B8
Target Name	Canine Coronavirus N Protein
Species	Virus
Immunogen	Full Virus of Canine Coronavirus
Conjugation	Un-conjugated
Alternate Names	COVID-19 virus, SARS-CoV-2 virus, coronavirus disease 2019 virus, severe acute respiratory syndrome coronavirus 2, 2019-nCoV.

Application Instructions

Application table	Application	Dilution
	ICC/IF	1:250 - 1:750
	WB	1:1000 - 1:1500
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	


Properties

Form	Liquid
Purification	Affinity purified.
Buffer	PBS, 0.05% Sodium azide and 20% Glycerol.
Preservative	0.05% Sodium azide
Stabilizer	20% Glycerol
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.



ARG67203 anti-Canine Coronavirus N Protein antibody [2B8] ICC/IF image

Immunofluorescence: ARG67203 anti-Canine Coronavirus N Protein antibody [2B8] were used for detecting Canine Coronavirus N Protein.



ARG67203 anti-Canine Coronavirus N Protein antibody [2B8] WB image

Western blot: ARG67203 anti-Canine Coronavirus N Protein antibody [2B8] were used for detecting Canine Coronavirus N Protein.