

ARG70026
anti-PRRSV antibody [SR30] (FITC)Package: 1 ml
Store at: -20°C

Summary

Product Description	FITC-conjugated Mouse Monoclonal antibody [SR30] recognizes PRRSV
Tested Reactivity	Virus
Tested Application	ICC/IF
Host	Mouse
Clonality	Monoclonal
Clone	SR30
Isotype	IgG1
Target Name	PRRSV
Species	Virus
Immunogen	A highly conserved epitope on the nucleocapsid protein of PRRSV.
Conjugation	FITC

Application Instructions

Application table	Application	Dilution
	ICC/IF	1:50
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

Properties

Form	Liquid
Buffer	PBS
Storage instruction	For continuous use, store undiluted antibody at 2-8°C in dark for up to a week. For long-term storage, aliquot and store at -20°C or below in dark. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.