

ARG70040
Human IgE protein (native) (Tag Free)

Package: Quote


Store at: -20°C

Summary

Product Description	Azide free Human monoclonal IgE, kappa light chain
Tested Application	SDS-PAGE
Host	Human
Clonality	Monoclonal
Clone	HE1
Isotype	IgE, kappa
Target Name	Human IgE
Species	Human
Conjugation	Un-conjugated

Properties

Form	Liquid
Purification	Purification with Protein L.
Purification Note	Endotoxin level is less than 25 EU/mg.
Purity	> 90% (by SDS-PAGE)
Buffer	PBS (pH 7.4)
Concentration	1 mg/ml
Storage instruction	For continuous use, store undiluted protein at 2-8°C for up to a week. For long term storage, aliquot and store at -20°C or -80°C. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening.
Note	For laboratory research only, not for drug, diagnostic or other use.



ARG70040 Human IgE protein (native) (Tag Free) SDS-PAGE image

SDS-PAGE analysis of ARG70040 Human IgE protein (native) (Tag Free) and IgE (purified from serum from a myeloma patient).