

ARG70219 SARS-CoV-2 Spike recombinant protein (RBD) (His-tagged, C-ter)

Package: 100 µg
Store at: -20°C


Summary

Product Description	E. coli expressed, His-tagged (C-ter) SARS-CoV-2 Spike recombinant protein (RBD)
Tested Application	SDS-PAGE
Target Name	SARS-CoV-2 Spike protein (RBD)
Species	Virus
A.A. Sequence	Arg319 - Asn541
Expression System	E. coli

Properties

Form	Powder
Purity	> 98% (by SDS-PAGE)
Buffer	PBS (pH 7.4)
Reconstitution	

Images



SARS-CoV-2 Spike recombinant protein (RBD)

ARG70219 SARS-CoV-2 Spike recombinant protein (RBD) (His-tagged, C-ter) SDS-PAGE image

SDS-PAGE analysis of ARG70219 SARS-CoV-2 Spike recombinant protein (RBD) (His-tagged, C-ter).