

ARG70221 HCoV-OC43 nucleocapsid recombinant protein (His-tagged, C-ter)

Package: 100 µg
Store at: -20°C


Summary

Product Description	E. coli expressed, His-tagged (C-ter) HCoV-OC43 nucleocapsid recombinant protein
Tested Application	SDS-PAGE
Target Name	HCoV-OC43 nucleocapsid protein
Species	Virus
A.A. Sequence	Met1 - Ile448
Expression System	E. coli

Properties

Form	Powder
Purity	> 98% (by SDS-PAGE)
Buffer	PBS (pH 7.4)
Reconstitution	

Images



ARG70221 HCoV-OC43 nucleocapsid recombinant protein (His-tagged, C-ter) SDS-PAGE image

SDS-PAGE analysis of ARG70221 HCoV-OC43 nucleocapsid recombinant protein (His-tagged, C-ter).