

## ARG70222 HCoV-HKU1 nucleocapsid recombinant protein (His-tagged, C-ter)

Package: 100 µg  
Store at: -20°C


### Summary

Product Description	E. coli expressed, His-tagged (C-ter) HCoV-HKU1 nucleocapsid recombinant protein
Tested Application	SDS-PAGE
Target Name	HCoV-HKU1 nucleocapsid protein
Species	Virus
A.A. Sequence	Met1 - Ala441
Expression System	E. coli

### Properties

Form	Powder
Purity	> 98% (by SDS-PAGE)
Buffer	PBS (pH 7.4)
Reconstitution	

### Images



HCoV-HKU1 nucleocapsid  
recombinant protein

ARG70222 HCoV-HKU1 nucleocapsid recombinant protein (His-tagged, C-ter) SDS-PAGE image

SDS-PAGE analysis of ARG70222 HCoV-HKU1 nucleocapsid recombinant protein (His-tagged, C-ter).