

ARG70356 SARS-CoV-2 ORF8 recombinant protein (His-SUMO-tagged, N-ter)

Package: 100 µg, 20 µg
Store at: -20°C


Summary

Product Description	HEK293 expressed, His-SUMO tagged (N-ter) SARS-CoV-2 ORF8 recombinant protein
Tested Application	SDS-PAGE
Target Name	SARS-CoV-2 ORF8 protein
Species	Virus
A.A. Sequence	Phe1 - Ile121
Expression System	HEK293

Properties

Form	Powder
Purity	> 95% (by SDS-PAGE)
Buffer	PBS (pH 7.4)
Reconstitution	

Images



SARS-CoV-2 ORF8
recombinant protein
(His-SUMO-tagged, N-ter)

ARG70356 SARS-CoV-2 ORF8 recombinant protein (His-SUMO-tagged, N-ter) SDS-PAGE image

SDS-PAGE analysis of ARG70356 SARS-CoV-2 ORF8 recombinant protein (His-SUMO-tagged, N-ter).