

## ARG70362 Monkeypox Virus L1R recombinant protein (His-tagged, C-ter)

Package: 100 µg  
Store at: -20°C


### Summary

Product Description	E. coli expressed, His-tagged (C-ter) Monkeypox Virus L1R recombinant protein
Tested Application	SDS-PAGE
Target Name	Monkeypox Virus L1R
Species	Virus
A.A. Sequence	Met1 - Gln152
Expression System	E. coli

### Properties

Form	Powder
Purity	> 95% (by SDS-PAGE)
Buffer	PBS (pH 7.4)
Reconstitution	

### Images



ARG70362 Monkeypox Virus L1R recombinant protein (His-tagged, C-ter) SDS-PAGE image

SDS-PAGE analysis of ARG70362 Monkeypox Virus L1R recombinant protein (His-tagged, C-ter).