

ARG70368

Monkeypox Virus M1R recombinant protein (His-SUMO tagged, N-ter)

Package: 100 µg

Store at: -20°C


Summary

Product Description	E. coli expressed, His-SUMO tagged (N-ter) Monkeypox Virus M1R recombinant protein
Tested Application	SDS-PAGE
Target Name	Monkeypox Virus M1R
Species	Virus
A.A. Sequence	Met1 - Asn250
Expression System	E. coli

Properties

Form	Powder
Purity	> 95% (by SDS-PAGE)
Buffer	PBS (pH 7.4)
Reconstitution	

Images



ARG70368 Monkeypox Virus M1R recombinant protein (His-SUMO tagged, N-ter) SDS-PAGE image

SDS-PAGE analysis of ARG70368 Monkeypox Virus M1R recombinant protein (His-SUMO tagged, N-ter)