

ARG70373 Influenza A H3N2 nucleocapsid recombinant protein (A/Hong Kong/1/1968) (His-tagged, N-ter)

Package: 100 µg
Store at: -20°C


Summary

Product Description	E. coli expressed, His-tagged (N-ter) Influenza A H3N2 nucleocapsid recombinant protein (A/Hong Kong/1/1968)
Tested Application	SDS-PAGE
Target Name	Influenza A H3N2 nucleocapsid protein
Species	Virus
A.A. Sequence	Met1 - Asn498
Expression System	E. coli

Properties

Form	Powder
Purity	> 95% (by SDS-PAGE)
Buffer	PBS (pH 7.4)
Reconstitution	

Images



ARG70373 Influenza A H3N2 nucleocapsid recombinant protein (A/Hong Kong/1/1968) (His-tagged, N-ter) SDS-PAGE image

SDS-PAGE analysis of ARG70373 Influenza A H3N2 nucleocapsid recombinant protein (A/Hong Kong/1/1968) (His-tagged, N-ter)