

ARG70374 Influenza A H1N1 nucleocapsid recombinant protein (A/New York/1682/2009) (His-tagged, N-ter)

Package: 100 µg
Store at: -20°C


Summary

Product Description	E. coli expressed, His-tagged (N-ter) Influenza A H1N1 nucleocapsid recombinant protein (A/New York/1682/2009)
Tested Application	SDS-PAGE
Target Name	Influenza A H1N1 nucleocapsid protein
Species	Virus
Expression System	E. coli

Properties

Form	Powder
Purity	> 95% (by SDS-PAGE)
Buffer	PBS (pH 7.4)
Reconstitution	

Images



ARG70374 Influenza A H1N1 nucleocapsid recombinant protein (A/New York/1682/2009) (His-tagged, N-ter) SDS-PAGE image

SDS-PAGE analysis of ARG70374 Influenza A H1N1 nucleocapsid recombinant protein (A/New York/1682/2009) (His-tagged, N-ter)