

ARG70376 Enterovirus A71 VP1 protein (His-tagged, C-ter)

Package: 100 µg
Store at: -20°C


Summary

Product Description	E. coli expressed, His-tagged (C-ter) Enterovirus 71 VP1 nucleocapsid recombinant protein.
Tested Application	SDS-PAGE
Target Name	Enterovirus 71
Species	Virus
A.A. Sequence	Gly1 - Ser297
Expression System	E. coli
Alternate Names	Enterovirus 71; Enterovirus A71; EV71; EV-A71; EV71 virus

Properties

Form	Powder
Purity	> 95% (by SDS-PAGE)
Buffer	PBS (pH 7.4)
Reconstitution	

Images



ARG70376 Enterovirus A71 VP1 protein (His-tagged, C-ter) SDS-PAGE image

SDS-PAGE analysis of ARG70376 Enterovirus A71 VP1 protein (His-tagged, C-ter).