

ARG70398 Mouse IL38 recombinant protein (His-tagged, C-ter)

Package: 100 µg, 20 µg
Store at: -20°C

Summary

Product Description	E. coli expressed, His-tagged (C-ter) Mouse IL38 recombinant protein
Tested Application	SDS-PAGE
Target Name	IL38
Species	Mouse
A.A. Sequence	Met1 - Arg152
Expression System	E. coli
Alternate Names	IL1F10; Interleukin 1 Family Member 10; IL-1HY2; IL1-Theta; IL-1F10; FKSG75; Interleukin-1 Receptor Antagonist FKSG75; Interleukin-1 Family Member 10; Family Of Interleukin 1-Theta; Interleukin-1 HY2; Interleukin-38; FIL1 Theta; IL-1 Theta; MGC119832; MGC119833; MGC11983; IL1HY2; IL-38; Interleukin-1 Receptor Antagonist-Like FIL1 Theta; Interleukin 1 Family, Member 10 (Theta); Interleukin 1 Family Member 10 (Theta); IL-1F10 (Canonical Form IL-1F10a); Interleukin-1 Theta; FIL1- Theta; FIL1-Theta; FIL1T; IL38


Properties

Form	Powder
Purification Note	Endotoxin level is less than 0.1 EU/µg of the protein, as determined by the LAL test.
Purity	> 98% (by SDS-PAGE)
Buffer	PBS (pH 7.4)
Reconstitution	It is recommended to reconstitute the lyophilized protein in sterile water to a concentration not less than 200 µg/mL and incubate the stock solution for at least 20 min at room temperature to make sure the protein is dissolved completely.
Storage instruction	For long term, lyophilized protein should be stored at -20°C or -80°C. After reconstitution, aliquot and store at -20°C or -80°C for up to one month. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Gene Symbol	IL1F10
Gene Full Name	Interleukin 1 Family Member 10
Background	The protein encoded by this gene is a member of the interleukin 1 cytokine family. This gene and eight other interleukin 1 family genes form a cytokine gene cluster on chromosome 2. This cytokine is thought to participate in a network of interleukin 1 family members to regulate adapted and innate immune responses. Two alternatively spliced transcript variants encoding the same protein have been reported.
Function	Cytokine with immunomodulatory activity. Alone, does not induce cytokine production, but reduces IL22 and IL17A production by T-cells in response to heat-killed Candida albicans. Reduces IL36G-induced production of IL8 by peripheral blood mononuclear cells. Increases IL6 production by dendritic cells stimulated by bacterial lipopolysaccharides (LPS). Ligand for IL-36R/IL1RL2.

Images



ARG70398 Mouse IL38 recombinant protein (His-tagged, C-ter) SDS-PAGE image

SDS-PAGE analysis of ARG70398 Mouse IL38 recombinant protein (His-tagged, C-ter)