

ARG80356 Human ANA Screen ELISA Kit

Package: 96 wells
Store at: 4°C

Summary

| | |
|---------------------|--|
| Product Description | ARG80356 Human ANA Screen ELISA Kit is an Enzyme Immunoassay kit for the determination of Human ANA Screen in serum and plasma (EDTA, heparin, citrate). |
| Tested Reactivity | Hu |
| Tested Application | ELISA |
| Target Name | ANA Screen |
| Conjugation | HRP |
| Sensitivity | diagn. 96.4 % |
| Sample Type | Serum and plasma (EDTA, heparin, citrate). |
| Standard Range | Cut - off |
| Sample Volume | 10 µl |

Application Instructions

| | |
|------------|------------------------|
| Assay Time | 30, 15, 15, 5 min (RT) |
|------------|------------------------|

Properties

| | |
|---------------------|--|
| Form | 96 well |
| Storage instruction | Store the kit at 2-8°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components. |
| Note | For laboratory research only, not for drug, diagnostic or other use. |

Bioinformation

| | |
|---------------|--|
| Highlight | Related products: ANA antibodies ; ANA ELISA Kits ; ANA recombinant proteins ; New ELISA data calculation tool: Simplify the ELISA analysis by GainData |
| Research Area | Immune System kit |