

## ARG80366 Human Ovarian antibody ELISA Kit

Package: 96 wells  
Store at: 4°C

### Summary

Product Description	ARG80366 Human Ovarian antibody ELISA Kit is a reliable and quantitative test for the determination of antibodies directed against human ovarian tissues. This test is intended for the use with serum.
Tested Reactivity	Hu
Tested Application	ELISA
Target Name	Ovarian antibody
Conjugation	HRP
Sensitivity	n.d
Sample Type	Serum.
Standard Range	6 - 100 U/ml
Sample Volume	5 µl

### Application Instructions

Assay Time	60, 60 min (37°C), 30 min (RT)
------------	--------------------------------

### Properties

Form	96 well
Storage instruction	Store the kit at 2-8°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.

### Bioinformation

Research Area	Cancer kit; Controls and Markers kit
---------------	--------------------------------------