

ARG80368 Human Anti-Spermatozoa antibody (ASA) ELISA Kit

Package: 96 wells

Store at: 4°C

Summary

Product Description	ARG80368 Human Anti-Spermatozoa antibody (ASA) ELISA Kit is a reliable and quantitative test for the determination antibodies directed against human spermatozoa. This test is intended for the use with seminal plasma. Please note: the terms "anti-spermatozoa antibodies", "anti-sperm antibodies" and "sperm antibodies" are equivalent. In these instructions the rather unwieldy but correct term "anti-spermatozoa antibodies" is used.
Tested Reactivity	Hu
Tested Application	ELISA
Target Name	Anti-Spermatozoa antibody (ASA)
Conjugation	HRP
Sensitivity	25 U/ml
Sample Type	Seminal plasma.
Standard Range	31 - 250 U/ml
Sample Volume	100 µl

Application Instructions

Assay Time	1 h, 1 h (37°C), 30 min (RT)
------------	------------------------------

Properties

Form	96 well
Storage instruction	Store the kit at 2-8°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Research Area	Immune System kit
---------------	-------------------