

ARG80372 Human C1q antibody ELISA Kit

Package: 96 wells
Store at: 4°C

Summary

Product Description	ARG80372 Human C1q antibody ELISA Kit is an Enzyme Immunoassay kit for the quantification of Human C1q Ab in serum and plasma (EDTA, heparin, citrate).
Tested Reactivity	Hu
Tested Application	ELISA
Target Name	C1q antibody
Conjugation	HRP
Sensitivity	0.5 U/ml
Sample Type	Serum and plasma (EDTA, heparin, citrate).
Standard Range	6.3 - 100 U/ml
Sample Volume	10 µl

Application Instructions

Assay Time	30, 15, 15, 5 min (RT)
------------	------------------------

Properties

Form	96 well
Storage instruction	Store the kit at 2-8°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Highlight	Related products: Complement component antibodies: Complement component ELISA Kits: Complement component Duos / Panels: New ELISA data calculation tool: Simplify the ELISA analysis by GainData
Research Area	Immune System kit; Neuroscience kit