

## ARG80390 Human GBM antibody ELISA Kit

Package: 96 wells  
Store at: 4°C

### Summary

Product Description	ARG80390 Human GBM antibody ELISA Kit is an Enzyme Immunoassay kit for the quantification of Human GBM Ab in serum and plasma (EDTA, heparin, citrate).
Tested Reactivity	Hu
Tested Application	ELISA
Target Name	GBM antibody
Conjugation	HRP
Sensitivity	1 U/ml
Sample Type	Serum and plasma (EDTA, heparin, citrate).
Standard Range	20 - 200 U/ml
Sample Volume	20 µl

### Application Instructions

Assay Time	30, 15, 15, 5 min (RT)
------------	------------------------

### Properties

Form	96 well
Storage instruction	Store the kit at 2-8°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.

### Bioinformation

Highlight	Related products:  New ELISA data calculation tool: <a href="#">Simplify the ELISA analysis by GainData</a>
Research Area	Cell Biology and Cellular Response kit; Developmental Biology kit