


ARG80483 Human Tryptophan ELISA Kit

Package: 96 wells
Store at: 4°C

Summary

| | |
|---------------------|--|
| Product Description | ARG80483 Human Tryptophan ELISA Kit is an Enzyme Immunoassay kit for the quantification of Tryptophan in serum, plasma (EDTA) and urine. |
| Tested Reactivity | Hu |
| Tested Application | ELISA |
| Specificity | |

Images



ARG80483 Human Tryptophan ELISA Kit standard curve image

ARG80483 Human Tryptophan ELISA Kit results of a typical standard run with optical density reading at 450 nm.