

ARG80653 Equine Estrone 3-Sulfate ELISA Kit

Package: 96 wells
Store at: 4°C

Summary

Product Description	ARG80653 Equine Estrone 3-Sulfate ELISA Kit is a competitive enzyme immunoassay for the quantification of Estrone 3-Sulfate in mare serum.
Tested Reactivity	Hrs
Tested Application	ELISA
Target Name	Estrone 3-Sulfate
Conjugation	HRP
Sensitivity	0.14 ng/ml
Sample Type	Serum.
Standard Range	5 - 1000 ng/ml
Sample Volume	20 µl

Application Instructions


Assay Time	1 h (shaking), 30 min (shaking), 30 min
------------	---

Properties

Form	96 well
Storage instruction	Store the kit at 2-8°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Highlight	Related products: Estrone ELISA Kits: New ELISA data calculation tool: Simplify the ELISA analysis by GainData
Research Area	Microbiology and Infectious Disease kit



ARG80653 Equine Estrone 3-Sulfate ELISA Kit standard curve image

ARG80653 Equine Estrone 3-Sulfate ELISA Kit results of a typical standard.