

ARG80657 Mouse / Rat Progesterone ELISA Kit

Package: 96 wells
Store at: 4°C

Summary

Product Description	ARG80657 Mouse / Rat Progesterone ELISA Kit is a competitive immunoassay for the quantification of progesterone in rat and mouse serum and plasma.
Tested Reactivity	Ms, Rat
Tested Application	ELISA
Target Name	Progesterone
Conjugation	HRP
Conjugation Note	Substrate: TMB and read at 450 nm
Sensitivity	0.04 ng/ml
Sample Type	Serum and plasma.
Standard Range	0.4 - 100 ng/ml
Sample Volume	25 µl

Application Instructions


Assay Time	60, 30 min
------------	------------

Properties

Form	96 well
Storage instruction	Store the kit at 2-8°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Highlight	Related products: Progesterone antibodies ; Progesterone ELISA Kits ; Progesterone Duos / Panels ; New ELISA data calculation tool: Simplify the ELISA analysis by GainData
Research Area	Cancer kit; Gene Regulation kit; Neuroscience kit; Signaling Transduction kit



ARG80657 Mouse/Rat Progesterone ELISA Kit standard curve image

ARG80657 Mouse/Rat Progesterone ELISA Kit results of a typical standard.