

ARG80755
Rat Estradiol ELISA KitPackage: 96 wells
Store at: 4°C

Summary


Product Description	ARG80755 Rat Estradiol ELISA Kit is an Enzyme Immunoassay kit for the quantification of Rat Estradiol in serum.
Tested Reactivity	Rat
Tested Application	ELISA
Target Name	Estradiol
Conjugation	HRP
Conjugation Note	TMB at 450 nm
Sensitivity	2.5 pg/ml
Sample Type	Serum.
Standard Range	5 - 1280 pg/ml
Sample Volume	75 µl

Properties

Form	96 well
Storage instruction	Store the kit at 2-8°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Highlight	Related products: Estradiol antibodies ; Estradiol ELISA Kits ; New ELISA data calculation tool: Simplify the ELISA analysis by GainData
Research Area	Neuroscience kit



ARG80755 Rat Estradiol ELISA Kit standard curve image

ARG80755 Rat Estradiol ELISA Kit results of a typical standard run with optical density reading at 450 nm.