

ARG81023 Human DGP / Deamidated gliadin proteins IgG antibody ELISA Kit

Package: 96 wells
Store at: 4°C

Summary

| | |
|---------------------|---|
| Product Description | ARG81023 Human DGP / Deamidated gliadin proteins IgG antibody ELISA Kit is an Enzyme Immunoassay kit for the quantification of Human's DGP / Deamidated gliadin proteins IgG antibody in serum and plasma (EDTA, heparin, citrate). |
| Tested Reactivity | Hu |
| Tested Application | ELISA |
| Target Name | DGP / Deamidated gliadin proteins IgG |
| Conjugation | HRP |
| Conjugation Note | TMB at 450 nm |
| Sensitivity | 1 U/ml |
| Sample Type | Serum and plasma (EDTA, heparin, citrate). |
| Standard Range | 6.3 - 100 U/ml |
| Sample Volume | 10 µl |

Properties

| | |
|---------------------|--|
| Form | 96 well |
| Storage instruction | Store the kit at 2-8°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components. |
| Note | For laboratory research only, not for drug, diagnostic or other use. |

Bioinformation

| | |
|---------------|---|
| Highlight | Related products: DGP antibodies ; DGP ELISA Kits ; New ELISA data calculation tool: Simplify the ELISA analysis by GainData |
| Research Area | Cell Biology and Cellular Response kit; Immune System kit |