

ARG81028 Human DHEA (free) ELISA Kit

Package: 96 wells
Store at: 4°C

Summary


| | |
|---------------------|--|
| Product Description | ARG81028 Human DHEA (free) ELISA Kit is an Enzyme Immunoassay kit for the quantification of Human DHEA (free) in saliva. |
| Tested Reactivity | Hu |
| Tested Application | ELISA |
| Target Name | DHEA |
| Conjugation | HRP |
| Conjugation Note | TMB at 450 nm |
| Sensitivity | 3.7 pg/ml |
| Sample Type | Saliva. |
| Standard Range | 10 - 2560 pg/ml |
| Sample Volume | 100 µl |

Properties

| | |
|---------------------|--|
| Form | 96 well |
| Storage instruction | Store the kit at 2-8°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components. |
| Note | For laboratory research only, not for drug, diagnostic or other use. |

Bioinformation

| | |
|---------------|---|
| Highlight | Related products: DHEA antibodies ; DHEA ELISA Kits ; New ELISA data calculation tool: Simplify the ELISA analysis by GainData |
| Research Area | Signaling Transduction kit |



ARG81028 Human DHEA (free) ELISA Kit standard curve image

ARG81028 Human DHEA (free) ELISA Kit results of a typical standard run with optical density reading at 450 nm.