

ARG81030 Aflatoxin (total) ELISA Kit

Package: 96 wells
Store at: 4°C

Summary


Product Description	ARG81030 Aflatoxin (total) ELISA Kit is an Enzyme Immunoassay kit for the quantification of Aflatoxin (total) in food (extraction, dilution, centrifugation).
Tested Reactivity	Other
Tested Application	ELISA
Target Name	Aflatoxin
Conjugation	HRP
Conjugation Note	TMB at 450 nm
Sensitivity	0.015 ng/ml
Sample Type	Food (extraction, dilution, centrifugation).
Standard Range	0.05 - 1.5 ng/ml
Sample Volume	100 µl

Properties

Form	96 well
Storage instruction	Store the kit at 2-8°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Highlight	Related products: Aflatoxin antibodies ; Aflatoxin ELISA Kits ; New ELISA data calculation tool: Simplify the ELISA analysis by GainData
Research Area	Microbiology and Infectious Disease kit



ARG81030 Aflatoxin (total) ELISA Kit standard curve image

ARG81030 Aflatoxin (total) ELISA Kit results of a typical standard run with optical density reading at 450 nm.