

ARG81201 6X His tagged Protein ELISA Kit

Package: 96 wells
Store at: 4°C

Summary

Product Description	ARG81201 6X His tagged Protein ELISA Kit is an Enzyme Immunoassay kit for the quantification of 6X His tagged proteins or residues in cell or tissue lysates
Tested Reactivity	Other
Tested Application	ELISA
Target Name	6X His tag
Conjugation	HRP
Conjugation Note	Read at 450 nm.
Sensitivity	1 ng/ml (Protein); 98 pM (Residue)
Sample Type	Cell or tissue lysates
Standard Range	0.98 - 4000 ng/ml (Protein); 0.098 - 400 nM (Residue)
Sample Volume	50 µl

Application Instructions


Assay Time	3.5 hours
------------	-----------

Properties

Form	96 well
Storage instruction	Store the kit at 2-8°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.


Bioinformation

Highlight	Related products: 6X His tag antibodies ; 6X His tag ELISA Kits ; 6X His tag Duos / Panels ; New ELISA data calculation tool: Simplify the ELISA analysis by GainData
-----------	--



ARG81201 6X His tagged Protein ELISA Kit standard curve image

ARG81201 6X His tagged Protein ELISA Kit results of a typical standard run with optical density reading at 450 nm. (His-tagged protein)



ARG81201 6X His tagged Protein ELISA Kit standard curve image

ARG81201 6X His tagged Protein ELISA Kit results of a typical standard run with optical density reading at 450 nm. (His-tag residue)