

ARG81227 Protein Carbonyl ELISA Kit

Package: 96 wells
Store at: 4°C, -20°C

Summary

Product Description	ARG81227 Protein Carbonyl ELISA Kit is an Enzyme Immunoassay kit for the quantification of Protein Carbonyl in purified protein, plasma, serum, cell lysate, or tissue homogenate.
Tested Reactivity	Other
Tested Application	ELISA
Target Name	Protein Carbonyl
Conjugation	HRP
Conjugation Note	Read at 450 nm.
Sensitivity	10 µg/ml
Sample Type	Purified protein, plasma, serum, cell lysate, or tissue homogenate
Standard Range	0.375 - 7.5 nmol/mg
Sample Volume	100 µl

Application Instructions


Assay Time	~ 7 hours
------------	-----------

Properties

Form	96 well
Storage instruction	Store components at 4°C or -20°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Highlight	Related products: New ELISA data calculation tool: Simplify the ELISA analysis by GainData
-----------	--



ARG81227 Protein Carbonyl ELISA Kit standard curve image

ARG81227 Protein Carbonyl ELISA Kit results of a typical standard run with optical density reading at 450 nm.