


ARG81365 Prostaglandin E2 ELISA Kit

Package: 96 wells
Store at: 4°C

Summary

Product Description	ARG81365 Prostaglandin E2 ELISA Kit is an Enzyme Immunoassay kit for the quantification of Prostaglandin E2 in saliva, serum, urine and cell culture supernatants.
Tested Reactivity	Other
Tested Application	ELISA
Specificity	

Images



ARG81365 Prostaglandin E2 ELISA Kit standard curve image

ARG81365 Prostaglandin E2 ELISA Kit results of a typical standard run with optical density reading at 405 nm.