

## ARG81392 Cortisol ELISA Kit

Package: 96 wells  
Store at: 4°C

### Summary

Product Description	ARG81392 Cortisol ELISA Kit is an Enzyme Immunoassay kit for the quantification of Cortisol in dried fecal samples, serum, plasma (EDTA and heparin), tissue culture media and urine.
Tested Reactivity	Other
Tested Application	ELISA
Target Name	Cortisol
Conjugation	HRP
Conjugation Note	Substrate: TMB and read at 450 nm.
Sensitivity	30.3 pg/ml
Sample Type	Dried fecal samples, serum, plasma (EDTA and heparin), tissue culture media and urine.
Standard Range	100 - 3200 pg/ml
Precision	Intra-assay: 8.8% Inter-assay: 8.13%

### Application Instructions


Assay Time	~ 1.5 hours
------------	-------------

### Properties

Form	96 well
Storage instruction	Store the kit at 2-8°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.

### Bioinformation

Highlight	Related products: <a href="#">Cortisol antibodies</a> ; <a href="#">Cortisol ELISA Kits</a> ; New ELISA data calculation tool: <a href="#">Simplify the ELISA analysis by GainData</a>
-----------	---



ARG81392 Cortisol ELISA Kit standard curve image

ARG81392 Cortisol ELISA Kit results of a typical standard run with optical density reading at 450 nm.