

ARG81629
HDL and LDL / VLDL Assay KitPackage: 100 tests
Store at: -20°C

Summary

Product Description	ARG81629 HDL and LDL / VLDL Assay Kit is a detection kit for the quantification of HDL and LDL / VLDL in serum.
Tested Reactivity	Other
Tested Application	FuncSt
Target Name	HDL and LDL / VLDL
Conjugation Note	Read at 340 nm.
Sensitivity	5 mg/dL
Sample Type	Serum.
Sample Volume	20 µl

Application Instructions

Application Note	Please note that this kit does not include a microplate.
Assay Time	~ 1 hour

Properties

Form	Liquid
Storage instruction	Store the kit at -20°C. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.