


ARG82001 Substance P ELISA Kit

Package: 96 wells
Store at: 4°C

Summary

Product Description	ARG82001 Substance P ELISA Kit is an Enzyme Immunoassay kit for the quantification of Substance P in serum and plasma samples.
Tested Reactivity	Hu, Ms, Rat
Tested Application	ELISA
Specificity	

Images



ARG82001 Substance P ELISA Kit standard curve image

ARG82001 Substance P ELISA Kit results of a typical standard run with optical density reading at 450 nm.