

## ARG82139 ATPase + GTPase Activity Assay Kit

Package: 200 tests  
Store at: 4°C

### Summary

Product Description	ARG82139 ATPase + GTPase Activity Assay Kit can be used to measure ATPase and GTPase activity.
Tested Reactivity	Hu, Ms, Rat, All
Tested Application	FuncSt
Target Name	ATPase + GTPase
Conjugation Note	Read at 620 nm
Sensitivity	0.007 U/L
Sample Volume	40 µl


### Application Instructions

Application Note	Please note that this kit does not include a microplate.
Assay Time	30 min

### Properties

Form	Liquid
Storage instruction	Store the kit at 2-8°C. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.

### Images



ARG82139 ATPase + GTPase Activity Assay Kit typical data demonstration image

ARG82139 ATPase + GTPase Activity Assay Kit results of a typical data with optical density reading at 620 nm.