

ARG82141
Bilirubin Assay Kit (Colorimetric)Package: 180 tests
Store at: 4°C

Summary


| | |
|---------------------|--|
| Product Description | ARG82141 Bilirubin Assay Kit (Colorimetric) is a detection kit for the quantification of Bilirubin in serum. |
| Tested Reactivity | Hu, Ms, Rat, All |
| Tested Application | FuncSt |
| Target Name | Bilirubin |
| Conjugation Note | Read at 530 nm |
| Sensitivity | 0.16 mg/dL |
| Detection Range | 1 - 10 mg/dL |
| Sample Type | serum |
| Sample Volume | 50 µl |

Application Instructions

| | |
|------------------|--|
| Application Note | Please note that this kit does not include a microplate. |
| Assay Time | 10 min |

Properties

| | |
|---------------------|---|
| Form | Liquid |
| Storage instruction | Store the kit at 2-8°C. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components. |
| Note | For laboratory research only, not for drug, diagnostic or other use. |



ARG82141 Bilirubin Assay Kit (Colorimetric) typical data demonstration image

ARG82141 Bilirubin Assay Kit (Colorimetric) results of a typical data with optical density reading at 530 nm.