

ARG82143 Chloride Assay Kit (Colorimetric)

Package: 250 tests
Store at: 4°C

Summary


Product Description	ARG82143 Chloride Assay Kit (Colorimetric) is a detection kit for the quantification of Chloride in biological, food and environmental samples.
Tested Reactivity	Hu, Ms, Rat, All
Tested Application	FuncSt
Target Name	Chloride
Conjugation Note	Read at 610 nm
Sensitivity	0.7 mg/dL (0.2 mM)
Detection Range	0.7 mg/dL (0.2 mM) - 35 mg/dL (10 mM)
Sample Type	biological, food and environmental samples
Sample Volume	5 µl

Application Instructions

Application Note	Please note that this kit does not include a microplate.
Assay Time	5 min

Properties

Form	Liquid
Storage instruction	Store the kit at 2-8°C. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.



ARG82143 Chloride Assay Kit (Colorimetric) typical data demonstration image

ARG82143 Chloride Assay Kit (Colorimetric) results of a typical data with optical density reading at 610 nm.