

**ARG82145**  
**Chromium Assay Kit (Colorimetric)**Package: 250 tests  
Store at: 4°C, -20°C

### Summary


Product Description	ARG82145 Chromium Assay Kit (Colorimetric) is a detection kit for the quantification of Chromium in serum, plasma, environmental, food and beverage.
Tested Reactivity	Hu, Ms, Rat, All
Tested Application	FuncSt
Target Name	Chromium
Conjugation Note	Read at 480 nm
Sensitivity	20 µg/L (0.38 µM, 0.02 ppb)
Detection Range	20 - 2000 µg/L
Sample Type	serum, plasma, environmental, food and beverage
Sample Volume	250 µl

### Application Instructions

Application Note	Please note that this kit does not include a microplate.
Assay Time	20 min

### Properties

Form	Liquid
Storage instruction	Store components at 4°C or -20°C. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.



ARG82145 Chromium Assay Kit (Colorimetric) typical data demonstration image

ARG82145 Chromium Assay Kit (Colorimetric) results of a typical data with optical density reading at 480 nm.