

**ARG82155**  
**Free Amino Nitrogen Assay Kit (Colorimetric)**Package: 100 tests  
Store at: 4°C

### Summary


Product Description	ARG82155 Free Amino Nitrogen Assay Kit (Colorimetric) is a detection kit for the quantification of Free Amino Nitrogen in food and beverages.
Tested Reactivity	Hu, Ms, Rat, All
Tested Application	FuncSt
Target Name	Free Amino Nitrogen
Conjugation Note	Read at 575 nm
Detection Range	0.2 - 10 mM
Sample Type	food and beverages
Sample Volume	5 µl

### Application Instructions

Application Note	Please note that this kit does not include a microplate.
Assay Time	20 min

### Properties

Form	Liquid
Storage instruction	Store the kit at 2-8°C. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.



ARG82155 Free Amino Nitrogen Assay Kit (Colorimetric) typical data demonstration image

ARG82155 Free Amino Nitrogen Assay Kit (Colorimetric) results of a typical data with optical density reading at 575 nm.