

ARG82160
Galactose Assay Kit (Colorimetric)

Package: 100 tests

Store at: -20°C

Summary


Product Description	ARG82160 Galactose Assay Kit (Colorimetric) is a detection kit for the quantification of Galactose in serum, plasma, urine, saliva, milk, culture medium, food and beverage.
Tested Reactivity	Hu, Ms, Rat, All
Tested Application	FuncSt
Target Name	Galactose
Conjugation Note	Read at 570 nm (colorimetric), or 530/585 nm (fluorometric)
Sensitivity	10 µM
Detection Range	10 - 1000 µM (colorimetric); 10 - 100 µM (fluorometric)
Sample Type	serum, plasma, urine, saliva, milk, culture medium, food and beverage
Sample Volume	20 µl

Application Instructions

Application Note	Please note that this kit does not include a microplate.
Assay Time	20 min

Properties

Form	Liquid
Storage instruction	Store the kit at -20°C. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.



ARG82160 Galactose Assay Kit (Colorimetric) typical data demonstration image

ARG82160 Galactose Assay Kit (Colorimetric) results of a typical data with optical density reading at 570 nm.