

ARG82164 Glucose 6 phosphate Assay Kit (Colorimetric)

Package: 100 tests

Store at: -20°C

Summary


Product Description	ARG82164 Glucose 6 phosphate Assay Kit (Colorimetric) is a detection kit for the quantification of Glucose 6 phosphate in serum, plasma, tissue and cell culture supernatants.
Tested Reactivity	Hu, Ms, Rat, All
Tested Application	FuncSt
Target Name	Glucose 6 phosphate
Conjugation Note	Read at 460 nm
Sensitivity	10 µM
Detection Range	10 - 1000 µM
Sample Type	serum, plasma, tissue and cell culture supernatants
Sample Volume	20 µl

Application Instructions

Application Note	Please note that this kit does not include a microplate.
Assay Time	20 min

Properties

Form	Liquid
Storage instruction	Store the kit at -20°C. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.



ARG82164 Glucose 6 phosphate Assay Kit (Colorimetric) typical data demonstration image

ARG82164 Glucose 6 phosphate Assay Kit (Colorimetric) results of a typical data with optical density reading at 460 nm.