

**ARG82168**  
**Glycolysis Assay Kit (Colorimetric)**

Package: 100 tests

Store at: -20°C

### Summary


Product Description	ARG82168 Glycolysis Assay Kit (Colorimetric) is a detection kit for the quantification of Glycolysis in cell culture supernatants.
Tested Reactivity	Hu, Ms, Rat, All
Tested Application	FuncSt
Target Name	Glycolysis
Conjugation Note	Read at 565 nm
Detection Range	3 - 10 mM
Sample Type	cell culture supernatants
Sample Volume	5 µl

### Application Instructions

Application Note	Please note that this kit does not include a microplate.
Assay Time	30 min

### Properties

Form	Liquid
Storage instruction	Store the kit at -20°C. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.



ARG82168 Glycolysis Assay Kit (Colorimetric) typical data demonstration image

Cell Density and Rate of Glycolysis: HL-60 cells seeded at varying cell densities in RPMI medium.