

## ARG82171 Invertase Activity Assay Kit (Colorimetric)

Package: 100 tests  
Store at: -20°C

### Summary


Product Description	ARG82171 Invertase Activity Assay Kit (Colorimetric) can be used to measure Invertase activity in biological and environmental samples.
Tested Reactivity	Hu, Ms, Rat, All
Tested Application	FuncSt
Target Name	Invertase
Conjugation Note	Read at 570 nm (colorimetric), or 530/585 nm (fluorometric)
Sensitivity	0.007 U/L
Sample Type	biological and environmental samples
Sample Volume	40 µl

### Application Instructions

Application Note	Please note that this kit does not include a microplate.
Assay Time	40 min

### Properties

Form	Liquid
Storage instruction	Store the kit at -20°C. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.



ARG82171 Invertase Activity Assay Kit (Colorimetric) typical data demonstration image

ARG82171 Invertase Activity Assay Kit (Colorimetric) results of a typical data with optical density reading at 570 nm.