

**ARG82176**  
**Lipase Activity Assay Kit (Colorimetric)**Package: 100 tests  
Store at: 4°C

### Summary


Product Description	ARG82176 Lipase Activity Assay Kit (Colorimetric) can be used to measure Lipase activity in Serum, plasma, saliva and urine.
Tested Reactivity	Hu, Ms, Rat, All
Tested Application	FuncSt
Target Name	Lipase
Conjugation Note	Read at 412 nm
Sensitivity	40 U/L
Detection Range	40 - 1600 U/L
Sample Type	Serum, plasma, saliva and urine
Sample Volume	10 µl

### Application Instructions

Application Note	Please note that this kit does not include a microplate.
Assay Time	20 min

### Properties

Form	Liquid
Storage instruction	Store the kit at 2-8°C. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.



ARG82176 Lipase Activity Assay Kit (Colorimetric) typical data demonstration image

ARG82176 Lipase Activity Assay Kit (Colorimetric) results of a typical data with optical density reading at 412 nm.