

## ARG82185 Nitric Oxide Synthase Inhibitor Screening Kit

Package: 100 tests  
Store at: -20°C

### Summary

Product Description	ARG82185 Nitric Oxide Synthase Inhibitor Screening Kit is a screening kit for inhibitor screening and evaluation of nitric oxide synthase inhibitors.
Tested Reactivity	Hu, Ms, Rat, All
Tested Application	FuncSt
Target Name	Nitric Oxide Synthase Inhibitor
Conjugation Note	Read at 540 nm
Sample Volume	10 µl


### Application Instructions

Application Note	Please note that this kit does not include a microplate.
Assay Time	2 - 3 hours

### Properties

Form	Liquid
Storage instruction	Store the kit at -20°C. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.

### Images



ARG82185 Nitric Oxide Synthase Inhibitor Screening Kit inhibitor titration curve image

Inhibitor titrations: iNOS from mouse was incubated with various concentrations of DPI. Each concentration of inhibitor contained 20 v/v% DMSO (final 2 v/v% in 50 µl reaction).