

## ARG82191 Phosphatase Activity Assay Kit (Colorimetric)

Package: 500 tests  
Store at: 4°C, -20°C

### Summary

Product Description	ARG82191 Phosphatase Activity Assay Kit (Colorimetric) can be used to measure Phosphatase activity.
Tested Reactivity	Hu, Ms, Rat, All
Tested Application	FuncSt
Target Name	Phosphatase
Conjugation Note	Read at 405 nm
Sensitivity	3 ng
Sample Volume	50 µl


### Application Instructions

Application Note	Please note that this kit does not include a microplate.
Assay Time	30 min

### Properties

Form	Liquid
Storage instruction	Store components at 4°C or -20°C. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.

### Images



ARG82191 Phosphatase Activity Assay Kit (Colorimetric) enzyme kinetics graph

Kinetics of 0.07 and 1.8 ng/µl Phosphatase reaction, using ARG82191 Phosphatase Activity Assay Kit (Colorimetric).