


ARG82210 cAMP ELISA Kit (Colorimetric)

Package: 96 wells
Store at: 4°C, -20°C

Summary


Product Description	ARG82210 cAMP ELISA Kit (Colorimetric) is an Enzyme Immunoassay kit for the quantification of cAMP in serum, plasma, urine, saliva, cell/tissue lysate and cell culture supernatants.
Tested Reactivity	Other
Tested Application	ELISA
Specificity	

Images



ARG82210 cAMP ELISA Kit (Colorimetric) standard curve image

ARG82210 cAMP ELISA Kit (Colorimetric) results of a typical standard run with optical density reading at 450 nm. (Non-Acetylated cAMP)



ARG82210 cAMP ELISA Kit (Colorimetric) standard curve image

ARG82210 cAMP ELISA Kit (Colorimetric) results of a typical standard run with optical density reading at 450 nm. (Acetylated cAMP)