


## ARG82211 cAMP ELISA Kit (Chemiluminescent)

Package: 96 wells  
Store at: 4°C, -20°C

### Summary


Product Description	ARG82211 cAMP ELISA Kit (Chemiluminescent) is an Enzyme Immunoassay kit for the quantification of cAMP in serum, plasma, urine, saliva, cell/tissue lysate and cell culture supernatants.
Tested Reactivity	Other
Tested Application	ELISA
Specificity	

### Images



ARG82211 cAMP ELISA Kit (Chemiluminescent) standard curve image

ARG82211 cAMP ELISA Kit (Chemiluminescent) results of a typical standard. (Non-Acetylated cAMP)



ARG82211 cAMP ELISA Kit (Chemiluminescent) standard curve image

ARG82211 cAMP ELISA Kit (Chemiluminescent) results of a typical standard. (Acetylated cAMP)