

**ARG82283**  
**NAD / NADH Assay Kit**Package: 96 wells  
Store at: 4°C, -20°C

### Summary

Product Description	ARG82283 NAD / NADH Assay Kit is a detection kit for the quantification of NAD / NADH content in serum, plasma, urine, tissue extracts, cell lysate, cell culture supernatants.
Tested Reactivity	Other
Tested Application	FuncSt
Target Name	NAD / NADH
Conjugation Note	Read at 570 nm.
Sensitivity	0.1 µmol/l
Detection Range	0.15 - 40 µmol/l
Sample Type	Serum, plasma, urine, tissue extracts, cell lysate, cell culture supernatants.
Sample Volume	20 µl

### Application Instructions

Assay Time	~ 10 min
------------	----------

### Properties

Form	96 well
Storage instruction	Store components at 4°C or -20°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.